



[Home](#) [Products](#) [What's New](#) [Ordering](#) [Technical Resources](#) [Contact Us](#) [About EPICENTRE](#)

## | LIGASES |

### T4 DNA Ligase, Cloned


[Order](#) | [Protocol](#) | [Related Products](#) | [Previous Item/Next Item](#)

T4 DNA Ligase is the most versatile and commonly used ligase for DNA cloning. This ATP-dependent enzyme covalently joins 5'-phosphates to 3'-hydroxylated termini at the blunt or compatible cohesive ends of double-stranded DNA fragments produced by restriction enzyme digestion or other enzymatic processes.<sup>1,2</sup> T4 DNA Ligase has no activity on single-stranded nucleic acids. Following a ligation reaction, T4 DNA Ligase may be inactivated by incubation at 65°C for 10 minutes.

EPICENTRE offers T4 DNA Ligase at both standard (2 U/μl) and high (10 U/μl) concentrations. High concentration T4 DNA Ligase is useful for obtaining maximum efficiency in blunt-end ligations. T4 DNA Ligase is supplied with a 10X Reaction Buffer and a 25 mM ATP Solution.

### Applications

- Ligation of blunt or cohesive-ended DNA fragments ([Figure 1](#)).
- Repair of nicks in double-stranded nucleic acids.<sup>3</sup>



**Figure 1. Ligation activity of EPICENTRE's T4 DNA Ligase.** T4 DNA Ligase at 10 U/μl was diluted sequentially in 10-fold increments. One microliter of each dilution of T4 DNA Ligase was incubated with Hind III-cut lambda DNA at 16°C in 1X Reaction Buffer containing ATP. Lane 1, kb ladder; Lane 2, no enzyme; Lane 3, 1 U/μl; Lane 4, 0.1 U/μl; Lane 5, 0.01 U/μl; Lane 6, 0.001 U/μl; Lane 7, 0.0001 U/μl; Lane 8, 0.00001 U/μl. Note that full enzyme activity is seen down to one one-thousandth unit; enzyme activity is still apparent even at one ten-thousandth unit.

**Unit Definition:** One Weiss unit of T4 DNA Ligase converts 1 nmole of <sup>32</sup>P from pyrophosphate into Norit-adsorbable material in 2 minutes at 37°C in 33 mM Tris-acetate (pH 7.8), 66 mM potassium acetate, 10 mM magnesium acetate, 0.5 mM DTT, and 1 mM ATP.<sup>4</sup> One Weiss unit equals approximately 100 cohesive-end ligation units.

**Storage Buffer:** 50% glycerol containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl, 0.1 mM EDTA, 0.1% Triton X-100, and 1 mM DTT

**T4 DNA Ligase 10X Reaction Buffer:** 330 mM Tris-acetate (pH 7.8), 660 mM potassium acetate, 100 mM magnesium acetate, and 5 mM DTT. The Reaction Buffer does not contain ATP, which must be added to the reaction to a final concentration of 0.5 - 1.0 mM. A 25 mM solution of ATP is included.

**Quality Control:** T4 DNA Ligase is functionally tested in cloning assays and is free of detectable contaminating DNA exo- and endonuclease and RNase activities.

### References

1. Helfman, D.M. *et al.* (1987) *Meth. Enzymol.* **152**, 349.
2. Wu, R. *et al.* (1987) *Meth. Enzymol.* **152**, 343.
3. Sambrook, J. *et al.* (1989) in: *Molecular Cloning: A Laboratory Manual* (2nd ed.), Cold Spring Harbor Laboratory Press, New York.
4. Weiss, B. *et al.* (1968) *J. Biol. Chem.* **243**, 4543.

To order EPICENTRE products, choose your country of origin. If your country is not listed, [Contact Us](#).

Australia

Catalog No.	Concentration	Size
<b>T4 DNA Ligase, Cloned</b>		
L0805H	2 U/μl	500 U
L0810H	2 U/μl	1,000 U
L0820H	2 U/μl	2,000 U
LH805H	10 U/μl	500 U
LH810H	10 U/μl	1,000 U
LH820H	10 U/μl	2,000 U

Includes 10X Reaction Buffer and a separate 25 mM ATP Solution.

T4 DNA Ligase is also available in bulk. Please inquire.

**You may wish to consider the following related products:**

[Ampligase® Thermostable DNA Ligase](#)

[Colony Fast-Screen™ Kit \(PCR Screen\)](#) **NEW**

[Colony Fast-Screen™ Kit \(Size Screen\)](#)

[Fast-Link™ DNA Ligation and Screening Kit](#)

[Fast-Link™ DNA Ligation Kit](#)

[HK™ Thermolabile Phosphatase](#)

[Mung Bean Nuclease](#)

[T4 DNA Polymerase](#)

[T4 Polynucleotide Kinase, Cloned](#)

[T4 RNA Ligase](#)

[\[top of page\]](#)

[Home](#) | [Products](#) | [What's New](#) | [Ordering](#) | [Technical Resources](#) | [Contact Us](#) | [About EPICENTRE](#)

**EPICENTRE®—Reagents and Kits for Genomics, Proteomics, and RNA Research**  
©2004 EPICENTRE Technologies. All Rights Reserved.

L Number	Hits	Search Text	DB	Time stamp
1	2	("6143527").PN.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:46
2	0	((("6143527").PN.) and random	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 12:45
3	0	((("6143527").PN.) and random\$2	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 12:45
4	2	((("6143527").PN.) and (fragment)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 12:54
5	2	((("6143527").PN.) and (fragment)) and DNase	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 13:31
6	2	6117679.pn.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:06
7	0	6117679.pn. and (flap)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 13:31
8	0	09840861.ap.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 14:52
9	4	840861.ap.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 14:53
10	1	840861.ap. and artificial	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 14:56
11	2	840861.ap. and restricted adj1 bank	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 14:57
12	1	6117679.pn. and DNaseI	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:06
13	1	(6117679.pn. and DNaseI) and (random\$4)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:06
14	1	((("6143527").PN.) and marker	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:24

15	1	((6117679.pn. and DNASEI) and (random\$4)) and marker	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:27
16	0	((6117679.pn. and DNASEI) and (random\$4)) and marker) and flag	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:40
17	1	6117679.pn. and ligat\$4	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:42
18	1	6117679.pn. and ligas\$4	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:43
19	1	6117679.pn. and liga\$6	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 15:58
20	2	840861.ap. and marker	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:01
21	110	(isolat\$4 adj1 (dna polynucleotide)) near5 (biotin marker heptin)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:15
22	11	stemmer.in. and flap	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:18
23	5085	overlap\$4 near2 flap	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:16
24	6	(overlap\$4 near2 flap) and (shuffl\$4)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:17
25	14	(overlap\$4 near2 flap) and (flap near5 cleav\$5)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:18
26	123	stemmer.in. and overlap\$5	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:19
27	92	(stemmer.in. and overlap\$5) and willem.in.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:20
28	90	((stemmer.in. and overlap\$5) and willem.in.) and (cleav\$5 digest\$5)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:21

29	4	((stemmer.in. and overlap\$5) and willem.in.) and ((cleav\$5 digest\$5) near5 (single adj1 strand (ss)))	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:41
30	119	flap adj1 endonuclease	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:41
31	105	(flap adj1 endonuclease) and clon\$4	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:43
32	105	((flap adj1 endonuclease) and clon\$4) and (cleav\$4 fragment\$5 digest\$5)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:44
33	54	((flap adj1 endonuclease) and clon\$4) and (flap near5 (cleav\$4 fragment\$5 digest\$5))	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:45
34	1	((("6143527").PN.) and (single adj1 strand\$5)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:49
35	2	6117679.pn. and (single adj1 strand\$4)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 16:49
36	1	((("6143527").PN.) and overlap\$4	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:25
37	7563	CDR and (liga\$5)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:24
38	7	(CDR and (liga\$5)) and (flap adj1 endonuclease)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:24
39	1	(flap adj1 endonuclease) near3 (cleav\$4 fragment\$ digest\$4) adj3 (nucleic)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:33
40	1	("20030215800").PN.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:34
41	1	((("20030215800").PN.) and flap	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:41
42	110	coco.in.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:41

43	1	coco.in. and levinson.in.	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:41
44	0	((("6143527").PN.) and (nick)	USPAT; US-PGPUB; EPO; JPO; DERWENT; IBM_TDB	2004/08/16 17:46